

Gymna FT - 100 Traction Device 161007

Datasheet



Cervical traction devices gently stretch the neck to reduce pressure on the spine by pulling or separating the vertebrae. The FT - 100 can be used for continuous and pulsed traction and comes with a patient emergency button.

- LCD touchscreen
- User definable force
- Programmable via USB
- All therapy parameters displayed on screen
- Compatible with gymna.GO traction table
- Adjustable screen contrast
- Adjustable volume

SPECIFICATIONS

Graphic LCD display:	320 x 240 pixels touchscreen	Power Supply:	230 VAC / 50-60 Hz / 55W
Operation modes:	Intermittent and Static	Standard accessories:	Patient safety button Power cable
High-Low Force Range:	1-100 kg / 1-99 kg \pm 2 kg	Device classification:	Ila
Speed Range:	8 levels (from 1:1 to 1:8)	Conformity:	IEC/EN 60601-1 IEC/EN 60601-1-2 MDD 93/42/EEC 2012/19/EC 2011/65/EU
Hold & Rest Time Range:	0 - 99 seconds		
Therapy Duration:	1 - 99 minutes		
Traction Rope Length:	160 cm		
Dimensions:	240 (W) x 340 (D) x 170 (H) mm		
Weight:	12 kg		