

Product Specification Sheet/ Gymna kontaktgel Medi-gel

General properties:

1.Appearance

Specification: transparant clear gel with no smell

Equipment: not applicable

Method: Examine the contents of at least one packaged product for homogeneity and the presence of foreign matter.

Evaluation: Note the colour as observed. There is no change permitted in uniformity and the presence of foreign matter is not acceptable.

2.pH

Specification: min 7.0 max 7.5

Equipment: pH meter

Method: Check that the pH meter is calibrated and when necessary calibrate in accordance with the operating instructions.

Evaluation: Note the reading on the pH meter.

3.Density

Specification: 1.0075-1.0125

Equipment: Densitometer

Method: Before using the densitometer read the operating instructions. Apply the sample as instructed. Carry out the measurement in duplicate.

Evaluation: Note the reading on the densitometer.

4.Microbiology:

Specification: Total plate count max. 500/g, pathogens absent.

Equipment: Analysis is carried out by a contract laboratory.

Method: Eur.Pharm. 4th ED.2002 Par.5.1.4.cat.2

5.Storage and transportation

In the original unopened containers the product can be stored for at least 30 months, protected from moisture at below 30 gr C.