


Art. 250174 D - Height adjustment

Tilt Table with two electric motors for height and tilt adjustment



- Height adjustment range 49 Δ 94 cm for easy transfer of wheelchair patients
- Gas spring assisted regulation of backrest
- The design of the support frame allows smooth lifting and transfer of the patient


- 2 electric motors with a thrust force **10.000 N** for height and **10.000 N** for tilt adjustment
- Electric height adjustment range 49 Δ 94 cm
- Independently adjustable and folding twin footboard. Positive/negative tilt + - 15°
- Tilt adjustment up to + 88°
- Fire-retardant padding and covering class 11M
- Gas spring assisted backrest tilt adjustment + 80°
- Equipped with set of **5** fixation belts
- Two fixed and **two** swivelling castors **with** brake
- Goniometer
- Vertical and horizontal safety work load: 160 kg
- Dimensions: 200 x 70 x 49 Δ 94 h cm



OPTIONS



250176 D



250175 D

B ACTUAL VERSION
RAFFAELLO

